

MONTHLY DEFIBRILLATOR CHECKLIST

To reorder accessories or AED contact First Aid Accident & Emergency **Ph: 07 55205068**

LOCATION OF AED		
MAKE/MODEL	EXPIRY OF PADS	
SERIAL NUMBER	EXPIRY OF BATTERIES	



CHECK

Check AED indicator for readiness

Check AED, wall cabinet or tough case for damage

Check expiry of battery

Check expiry of pads
(if not combined with battery)

Alarmed Wall Cabinet 9V battery:
Replace on 1 April or annually

Check AED prep pack for inclusions:
razor, shears, CPR face shield, gloves

Check all AED signage is visible

DATE CHECKED

INITIALS

[illegible]